

Use of WHOOP devices at ITF events

All WHOOP devices (3.0, 4.0, 5.0, MG) are approved for use at all ITF events, as per the Duties and Procedures (Interpretations III, Electronic Devices, Pg 41).

WHOOP 4.0, 5.0, and MG devices have haptic feedback functions that must be disabled during matches.

Haptic feedback is restricted during matches for all PAT devices to prevent you or your opponent receiving messages during points.

Players are responsible for ensuring these settings are correct before they walk onto court.

If players are found using their WHOOP device in a prohibited mode, the player can be penalised under the code of conduct for coaching and will be asked to remove their WHOOP.

Instructions for how to configure settings correctly are detailed in the linked WHOOP PAT Report.



Link to full WHOOP PAT Report

